



## SmartVision and SmartVision LITE


### Release note – software version n°1454.2664

#### Software Update Process

SmartVision regularly checks the availability of updates on your device and you will be notified automatically on home screen with a red label on the “Settings” icon. .

However, you can manually check the availability an update from "Settings"  menu, then select “About phone”, and finally select “System update”

To start the update process:

- select « Settings » item from home screen, 
- then select “Download” item
- Download will start automatically and will last for about 30 to 40 minutes depending of your network bandwidth. Subsequently, the device will reboot automatically until “Home screen” is displayed and announced.

#### Note:

- Please do not use SmartVision when the download is in progress.
- You must have data connection, Wi-Fi or 3G enabled to download the software update. Because of important size of this update (about 1GB), we recommend to use the Wi-Fi connection for downloading.
- If the device goes to sleep mode during the download process, just make a short press on ON/OFF button to wake-up the device, then when the download is complete, the upgrade process will start automatically.

#### SmartVision System improvements:

- **Overview:**
  - o Bug fixed and stability improvements.
  - o Translation fixed and improved in Italian, Spanish, German and Netherlands languages.
- **Accessibility:**
  - o Software upgrade process: download progression bar is now vocalized, and estimated “time remaining” information has been added.
  - o new settings have been added to confirm when a character entered by the keypad has been validated (see into Settings / Accessibility) :
    - Play a beep after the validation of a character.
    - Make a vibration after a validation of a character.
    - Validation timeout: 1, 2, 3, 4 or 5 seconds.

- **User Interface:**

- Speech speed setting value is now saved when using key shortcuts “#” + “4” or “#” + “5” (correction).
- Headset : Maximum output sound level has been increased for FM and Music applications.
- Edition mode : New keys attribution for a better use
  - Continuous press on star “★” key and short press on key “1”: moves the cursor one character to the left.
  - Continuous press on star “★” key and short press on key “2”: moves the cursor one character to the right.
  - Continuous press on star “★” key and short press on key “3”: copy the selected text.
  - Continuous press on star “★” key and short press on key “4”: moves the cursor one word to the left.
  - Continuous press on start “★” key and short press on key “5”: moves the cursor one word to the right.
  - Continuous press on star “★” key and short press on key “6”: cut the selected text.
  - Continuous press on star “★” key and short press on key “7”: moves the cursor back to the beginning of the text.
  - Continuous press on star “★” key and short press on key “8”: moves the cursor at the end of the text.
  - Continuous press on star “★” key and short press on key “9”: paste the selected text.
  - Continuous press on star “★” key and short press on pound “ # “ key : selects the whole text.

- **Vocalization:**

- ON/OFF screen (long press on ON/OFF key): title of the windows is now vocalized.
- Wrong format announcement for time and time duration information has been corrected.
- The name of the network operator and status of the data connection (3G, Edge...) are now announced when asking for the phone status (long press on key “#” followed by key “1”).

- **New features / applications :**

- ON / OFF screen: option to switch OFF permanently the display (black screen) in order to improve the battery life.
- Color detection application:
  - Color detection: point the camera on an object, and the color will be announced automatically. Press the OK key to repeat.
  - Find color: select a color from the list and when the camera will see this color, a sound is played.

## Applications:

- **Alarm:**
  - o The recurrence days are now fully vocalized (Monday, Tuesday ...)
- **Calculator**
  - o Correction of wrong vocalization of "Division" operator announced "percentage".
- **Camera**
  - o "Simple" display mode is set by default for an easier use.
- **Contact:**
  - o Issue to un-flag as favorite a contact imported from the SIM card is now corrected.
  - o A new option has been added on "Menu" to delete all contacts stored into the phone memory.
  - o Possibility to create a new contact directly from "Phone" and "Messaging" application.
- **Email:**
  - o Only contact including valid Email address are proposed as candidates when entering the destination name (correction).
  - o When writing a new Email, the "Send" button is now always displayed at the bottom of the screen.
  - o Multi-selection mode has been added in order to select and delete several E-mails at once.
  - o Improvement of the "text mode" by removing all E-mail's header for a better understanding.
- **FM Radio:**
  - o Warning pop-up has been added if the headset, used for the antenna, is not plugged on the device.
- **Gallery:**
  - o The name of the picture is now vocalized when selected.
- **Kapten (GPS):**
  - o Entry text area for city name is now vocalized when selected (correction)
  - o "You're arrived" message has been added at the end of the navigation when the destination has been reached, and at the end of the roadbook.
  - o A confirmation pop'up has been added before closing the application.
  - o On "Direct guidance" mode, navigation instructions are now repeated automatically (correction).
  - o key shortcut has been added during navigation:
    - Key "1": Where am I?
    - Key "2": Last instruction.
    - Key "3": Next instruction.
    - Key "4": Current speed.
    - Key "5": Speed limit.
    - Key "6": Remaining distance.
    - Key "7": Point Of Interest around.
    - Key "8": Enable or disable free navigation.
    - Key "9": GPS signal.

- **Messaging (SMS/MMS):**
  - Into a MMS, the saving of an attached item (picture, music, V-card) has been improved by using a long press on “OK” key on selected message.
- **OCR:**
  - Accessibility improvement to take a picture just by pressing “OK” key when “Manual mode” is selected.
  - Accessibility improvement to select a picture from gallery by using up and down key.
- **Phone:**
  - Into phone number text area, long press on “0” key insert the “+” character for international call (correction).
  - Output sound level management during a call has been improved.
  - During a call, the phone vibrates when the “hands free” mode is enabled or disabled by using a long press on voice command button.
- **SOS :**
  - A full “Red Warning” message is now displayed on the screen during SOS process.